

Nanosecond Laser Facility

Name of the facility:

Nanosecond Nd:YAG Pumped Dye Laser

Tunable Dye laser with output from 400 – 800 nm, with the option of frequency doubling. Complete tuning range 210 – 800 nm


Make and Model:

Radiant Dyes, Narrow Scan GR

Specification/ Features:

Tunable Dye laser with output from 400 – 800 nm, (with the option of frequency doubling)

Complete tuning range 210 – 800 nm.



Available mode of use: Based on individual requirements

Whether the facility is open to external users: YES

Location:

Room no: 206

Department: Chemistry

Contact details : 022-25767182 / 022-25764191

Convenor/in-charge name : G. Naresh Patwari

Email id : naresh@chem.iitb.ac.in